



ANJ BIOMEDICALS LLC Tel: 302.357.7443
630 Route 206, Bordentown,NJ , 08505
www.anjbiomedicals.com info@anjbiomedicals.com

CERTIFICATE OF ANALYSIS

Name Of Item : D-GALACTOSE

Formula : C₆H₁₂O₆

M.W. : 180.16

Cat. No. : 100122

CAS NO. : 59-23-4

Batch No. :

Date Of Mfg. :

Date of Analysis :

Type Of Test	Standard	Observed
Description	White powder	White powder
Purity	98 %	99%
Melting point	162 - 166°C	165.8°C
Specific rotation at 20°C	+79° TO +81°	+79.3°
Heavy metals (as Pb)	<0.001%	<0.001%
Glucose	<0.01%	<0.01%
Water	<0.05%	0%

This Document is electronically generated & hence valid without signature