



ANJ BIOMEDICALS LLC Tel: 302.357.7443
630 Route 206, Bordentown, NJ , 08505
www.anjbiomedicals.com
info@anjbiomedicals.com

CERTIFICATE OF ANALYSIS

PRODUCT NAME: CRYSTAL VIOLET (GENTION VIOLET)
PRODUCT NUMBER: SW3874
BATCH NUMBER: 546455
CAS NUMBER: 548-62-9
MOLECULAR FORMULA : C₂₅H₃₀N₃Cl
MOLECULAR WEIGHT : 407.98
MFG. DATE: OCT-2022
EXP. DATE: SEP-2024

| TEST | SPECIFICATIONS | RESULTS |
|---|----------------|------------|
| Appearance (Colour) | Dark green | Dark green |
| Appearance (Form) | Powder | Powder |
| Loss on drying | max. 10% | 5.5% |
| Absorbance (A) of 1% aq. solution in a 1cm cell v/s H ₂ O @590nm | 2000 - 2450 | 2370 |

This Document is electronically generated & hence valid without signature