



ANJ BIOMEDICALS LLC Tel: 302.357.7443
630 Route 206, Bordentown, NJ , 08505
www.anjbiomedicals.com
info@anjbiomedicals.com

CERTIFICATE OF ANALYSIS

PRODUCT NAME: D.P.X Mountant
PRODUCT NUMBER: 200656
BATCH NUMBER: 124179
CAS NO : 130-12-2
MFG. DATE: OCTOBER-2022
EXP. DATE: SEPTEMBER-2027

Type of Test	Specification	Results
Description	Clear colorless liquid	Clear colorless liquid
Color	Water white	Water white
Refractive index	1.515-1.525	1.520
Natural drying time	40 minutes	40 minutes

This Document is electronically generated & hence valid without signature