



ANJ BIOMEDICALS LLC Tel: 302.357.7443
630 Route 206, Bordentown, NJ, 08505
www.anjbiomedicals.com
info@anjbiomedicals.com

CERTIFICATE OF ANALYSIS

PRODUCT NAME: DIPHENYLAMINE AR
PRODUCT NUMBER: 200648
BATCH NUMBER: 962078
CAS NO: 122-39-4
M.W : 169.23
FORMULA : C₁₂H₁₁N
MFG. DATE: SEPT-2022
EXP. DATE: AUGUST-2027

Type of Test	Specification	Results
Description	Cream colored flakes/powder	Cream colored flakes/powder
Assay	Min 99%	99.77%
Min. sensitivity to nitrate (NO ₃)	1:1600000	1.1600000
Melting Point	53-55°C	53°C
Insoluble Matter	<0.01%	<0.01%
Sulphated ash	<0.02%	
Nitrate	No reaction	No reaction
Aniline(C ₆ H ₅ .NH ₂)	<0.01%	<0.01%

This Document is electronically generated & hence valid without signature