



**ANJ BIOMEDICALS LLC** Tel: 302.357.7443  
630 Route 206, Bordentown, NJ , 08505  
www.anjbiomedicals.com  
info@anjbiomedicals.com

## **CERTIFICATE OF ANALYSIS**

**PRODUCT NAME:** Sulfasalazine  
**PRODUCT NUMBER:** 200883  
**BATCH NUMBER:** 128996  
**CAS NO :** 599-79-1  
**M.W:** 398.39  
**FORMULA :** C<sub>18</sub>H<sub>14</sub>N<sub>4</sub>O<sub>5</sub>S  
**MFG. DATE:** OCTOBER-2022  
**EXP. DATE:** SEPTEMBER-2027

<b>Type of Test</b>	<b>Specification</b>	<b>Results</b>
Appearance (Form)	Powder	Powder
Appearance (Color) Yellow to brown	Conforms	Conforms
Identification A, B according to USP	Conforms	Conforms
Loss on Drying	≤1.0%	1.0%
Residue on Ignition	≤0.5%	0.5%
Chromatographic Purity(TLC) total impurity=4.0%	Conforms	Conforms
Assay(USP, on dried basis)	97.0-101.5%	101.5%

**This Document is electronically generated & hence valid without signature**