



ANJ BIOMEDICALS LLC Tel: 302.357.7443
630 Route 206, Bordentown, NJ , 08505
www.anjbiomedicals.com info@anjbiomedicals.com

CERTIFICATE OF ANALYSIS

PRODUCT NAME: HYDROXYPROPYL CELLULOSE
PRODUCT NUMBER: SW4498
BATCH NUMBER: C9848
CAS NUMBER: 9004-64-2
MFG. DATE: DEC.2022
EXP. DATE: NOV.2027

Test	Specification	Result
Appearance (Color)	White to Off-White	White
Appearance (Form)	Powder or Crystals	Powder
Infrared Spectrum	Conforms to Structure	Conforms
Viscosity c = 10%, Water, 25°C	250.0 - 800.0 cps	390.0 cps

This Document is electronically generated & hence valid without signature